

INTRAOCULAR LENSES



NOT AVAILABLE FOR SALE IN THE U.S.

VisiLens PMMA Single-Piece Posterior Chamber Intraocular Lens*

LENS DESIGN	MODEL NO.	OPTIC DIAMETER	OVERALL DIAMETER	LOOP	VAULT HEIGHT	A CONSTANT* ACD (mm)*	DIOPTER RANGE
	1E50C	5.0mm Equiconvex	12.0mm	Modified "C"	0.3mm	118.5 5.0	10 to 14 in 1.0 D 15 to 25 in 0.5 D 26 to 30 in 1.0 D
	1E55C	5.5mm Equiconvex	12.0mm	Modified "C"	0.3mm	118.5 5.0	
	1E60C	6.0mm Equiconvex 2 holes	12.5mm	Modified "C"	0.4mm	118.5 5.0	
	1E65C	6.5mm Equiconvex 2 holes	13.5mm	Modified "C"	0.4mm	118.5 5.0	

- Equiconvex UV-absorbing optics for less aberrations

-90° deviated haptics ensures self-centering and broad contact

-Step vaulted haptics for optimal stability

VisiLens PMMA Anterior Chamber Intraocular Lens**

LENS DESIGN	MODEL NO.	OPTIC DIAMETER	OVERALL DIAMETER	LOOP	A CONSTANT* ACD (mm)*	DIOPTER RANGE
	A1E283	5.5mm Equiconvex	12.8mm	Kelman Multiflex Design 3 point fixation 0.5mm step vault	115.3 3.0	10 to 14 in 1.0 D 15 to 25 in 0.5 D 26 to 30 in 1.0 D

** Custom-made lenses available upon request

VisiFold Hydrophilic Foldable Intraocular Lens

LENS DESIGN	MODEL NO.	OPTIC DIAMETER	OVERALL DIAMETER	LOOP	A CONSTANT* ACD (mm)*	DIOPTER RANGE
	F1C575	5.75mm	12.25mm	Broad "C"	118.0 5.0	10 to 14 in 1.0 D 15 to 25 in 0.5 D 26 to 30 in 1.0 D

*A Constant and ACD are estimated values only.

VisiRing Capular Tension Rings

RING DESIGN	MODEL NO.	OVERALL DIAMETER (mm)	CLOSED DIAMETER (mm)
	VisiRing 10	12.00	10.00
	VisiRing 11	13.00	11.00



VISIONARY MEDICAL SUPPLIES, INC.
SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE

6441 Enterprise Lane, Madison, WI 53719 USA

TEL: 608-270-3880

FAX: 608-270-3882

www.visionarymedicalsupplies.com