

OPHTHALMIC SUTURES

Visionary Medical sutures have undergone testing and received testimonials from ophthalmologists who found them to be of “equal or superior quality” to the sutures of leading competitors in every way. Our sutures are made from hardened steel and precisely formed to an exquisite mirror-like-tip with a square shank that is easy to hold and maneuver. Ophthalmologists find that the hard, sharp needle passes easily through tissue and resists deflection.



CE AND FDA

THREADS

Black Monofilament Nylon

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED	Advanced Micro Spatula Needle			
				CURVATURE	WIRE	LENGTH
10-0, 3/8 circle	N104062 (30cm)	N104061 (15cm)	ST140-6 3/8 circle (TG140-6)	140°	0.15mm	6.00mm
9-0, 3/8 circle	N94062 (30cm)	N94061 (15cm)				
10-0, 1/2 circle	N106062 (30cm)	N106061 (15cm)	ST160-6 1/2 circle (TG160-6)	160°	0.15mm	5.60mm
9-0, 1/2 circle	N96062 (30cm)	N96061 (15cm)				
10-0, 3/8 circle	N104082 (40cm)	N104081 (40cm)	ST140-8 3/8 circle (TG140-8)	140°	0.20mm	6.00mm
9-0, 3/8 circle	N94082 (40cm)	N94081 (40cm)				

NEEDLES

Black Braided Silk

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED	Advanced Micro Spatula Needle			
				CURVATURE	WIRE	LENGTH
8-0, 3/8 circle	S84082 (30cm)	S84081 (15cm)	ST140-8 3/8 circle (TG140-8)	140°	0.20mm	6.00mm

Virgin Silk (Not for use in the U.S.)

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED	Advanced Micro Spatula Needle			
				CURVATURE	WIRE	LENGTH
8-0, 3/8 circle	VS84082 (40cm)	VS84081 (40cm)	ST140-8 3/8 circle (TG140-8)	140°	0.20mm	6.00mm



VISIONARY MEDICAL SUPPLIES, INC.
SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE

6441 Enterprise Lane, Madison, WI 53719 USA

TEL: 608-270-3880

FAX: 608-270-3882

www.visionarymedicalsupplies.com