






OPHTHALMIC SUTURES

Visionary Medical sutures have undergone testing and received testimonials from ophthalmologists who found them to be “equal or superior quality” to the sutures of leading competitors in every way. Our sutures are made from hardened steel and precisely formed to an exquisite mirror-like-tip with a square shank that is easy to hold and maneuver. Ophthalmologists find that the hard, sharp needle passes easily through tissue and resists deflection.





ADVANCED MICRO SPATULA NEEDLE


Nylon, Black Monofilament

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
10-0, 3/8 circle	N104062 (30cm)	N104061 (15cm)	ST140-6 3/8 circle (TG140-6)	140°	0.15mm	6.00mm
9-0, 3/8 circle	N94062 (30cm)	N94061 (15cm)				
8-0, 3/8 circle	N84082 (30cm)*		ST140-8 3/8 circle (TG140-8)	140°	0.20mm	6.00mm
6-0, 3/8 circle	N640112 (76cm)*					
10-0, 1/2 circle	N106062 (30cm)	N106061 (15cm)	ST160-6 1/2 circle (TG160-6)	160°	0.15mm	5.60mm
9-0, 1/2 circle	N96062 (30cm)	N96061 (15cm)				
10-0, 3/8 circle	N104082 (40cm)	N104081 (40cm)	ST140-8 3/8 circle (TG140-8)	140°	0.20mm	6.00mm
9-0, 3/8 circle	N94082 (40cm)	N94081 (40cm)				


Silk, Black Braided

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
8-0, 3/8 circle	S84082 (30cm)	S84081 (15cm)	ST140-8 3/8 circle (TG140-8)	140°	0.20mm	6.00mm
6-0, 3/8 circle	S640112 (76cm)*					
6-0, 3/8 circle	S640112 (76cm)*		ST140-11 3/8 circle	140°	0.28mm	8.00mm
						

Silk, Virgin

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
8-0, 3/8 circle	VS84082 (40cm)	VS84081 (40cm)	ST140-8 3/8 circle (TG140-8)	140°	0.20mm	6.00mm
						

* FDA only



VISIONARY MEDICAL SUPPLIES, INC.

SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE

6441 Enterprise Lane, Madison, WI 53719 USA


TEL: 608-270-3880

FAX: 608-270-3882

www.visionarymedicalsupplies.com




OPHTHALMIC SUTURES

Visionary Medical sutures have undergone testing and received testimonials from ophthalmologists who found them to be "equal or superior quality" to the sutures of leading competitors in every way. Our sutures are made from hardened steel and precisely formed to an exquisite mirror-like-tip with a square shank that is easy to hold and maneuver. Ophthalmologists find that the hard, sharp needle passes easily through tissue and resists deflection.



ADVANCED MICRO SPATULA NEEDLE

Polypropylene, Black Monofilament

	MODEL-DOUBLE ARMED		CURVATURE	WIRE	LENGTH
9-0 straight	P90062 (20cm)*	ST000-6 straight 	0°	0.15mm	16.0mm
9-0, straight and 1/2 circle	P9006062 (30cm)*	ST000-6, and ST160-6 straight, 1/2 circle  	160°	0.15mm	16.0mm and 6.2mm


Polyglycolic Acid, Violet Braided

			CURVATURE	WIRE	LENGTH
6-0, 1/4 circle		ST90-1308 1/4 circle (NW2670) 	90°	0.33mm	8.0mm




TAPER POINT


Nylon, Black Monofilament

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
7-0, 3/8 circle	N74092 (30cm)*		TP140-9 3/8 circle (BV130-5) 	140°	0.23mm	6.00mm

Silk, Black Braided

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
7-0, 3/8 circle	S74092 (45cm)*		TP140-9 3/8 circle (BV130-5) 	140°	0.23mm	6.00mm

* FDA only



VISIONARY MEDICAL SUPPLIES, INC.

SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE

6441 Enterprise Lane, Madison, WI 53719 USA

TEL: 608-270-3880

FAX: 608-270-3882

www.visionarymedicalsups.com