

# VISI FOLD

## Hydrophilic Foldable Intraocular Lenses


VisiFold lenses are sold with the VisiGlide injector.  
Use VisiGlide to minimize surgical trauma and immediate postoperative inflammation.



NOT AVAILABLE FOR SALE IN THE U.S.



### DESCRIPTION



Material:	poly-HEMA
Refractive Index:	1.463
Optic Diameter:	5.75mm
Overall Diameter:	12.25mm
Water Content:	25%
Position Holes:	0
Angulation:	0 Degrees
Loops:	Broad "C"
'A' Constant:*	118.0
ACD:*	5.0

### POWER RANGE

10 to 15 in 1.0 Diopter
15 to 25 in 0.5 Diopter
25 to 30 in 1.0 Diopter

### CHARACTERISTICS

- Zero degree angulation facilitates error-free implantation
- Requires 2.8mm incision using VisiGlide
- Square edges prevent PCO
- Uniform center thickness across the entire power range
- Equiconvex UV absorbing optics for less aberrations
- Clinically proven safety and efficacy
- Specifications comply with ISO 11979 series of standards

\* 'A' Constant and ACD are estimated values only



**VISIONARY MEDICAL SUPPLIES, INC.**  
SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE

6441 Enterprise Lane, Madison, WI 53719 USA

TEL: 608-270-3880

FAX: 608-270-3882

[www.visionarymedicalsupplies.com](http://www.visionarymedicalsupplies.com)