


VISI RING

Capsular Tension Rings




NOT AVAILABLE FOR
SALE IN THE U.S.

VisiRing I0



OVERALL DIAMETER: 12mm
CLOSED DIAMETER: 10mm

VisiRing I1



OVERALL DIAMETER: 13mm
CLOSED DIAMETER: 11mm



ADVANTAGES

- Ensures stable conditions during surgery
- Circular expansion of capsular bag
- Prevents Capsular Fibrosis
- Better IOL centration when zonules are missing or damaged

INDICATIONS

- OL subluxation
- During usage of soft IOLs
- Damaged zonules, e.g. cases of Pseudoexfoliation or different degrees of wearing due to old age

PRECAUTIONS

- The leading end of the capsular tension ring may entangle with the bag fornix. This is indicated by traction folds in the posterior capsule. If the end is forcibly pushed farther, it may perforate the capsule.
- Entanglement occurs more readily with large or rigid capsular tension rings. Such entanglement can be avoided with adequate capsulorhexis sizing and with the use of high viscosity viscoelastic substances.



VISIONARY MEDICAL SUPPLIES, INC.

SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE

6441 Enterprise Lane, Madison, WI 53719 USA

TEL: 608-270-3880

FAX: 608-270-3882

www.visionarymedicalsupplies.com