


OPHTHALMIC SUTURE SPECIFICATIONS

ABOUT OUR SUTURES

Visionary Medical sutures are made from hardened steel and precisely formed to exquisite tips with a square shank that is easy to hold and maneuver. Ophthalmologists find that the hard, sharp needle passes easily through tissue and resists deflection. Our sutures have received testimonials from ophthalmologists who found them to be of "equal or superior quality" to the sutures of leading competitors in every way.



CE AND
FDA 510(k)



VISIONARY MEDICAL SUPPLIES, INC.









SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE

OPHTHALMIC SUTURE


SPECIFICATIONS

NYLON, BLACK MONOFILAMENT

Advanced Micro Spatula Needle

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED			CURVATURE	WIRE	LENGTH
10-0, 1/2 circle	N106062 (30cm)	N106061 (15cm)	ST160-6 1/2 circle (TG160-6)		160°	0.15mm	5.60mm
10-0, 3/8 circle	N104062 (30cm)	N104061 (15cm)	ST140-6 3/8 circle (TG140-6)		140°	0.15mm	6.00mm
10-0, 3/8 circle	N104082 (40cm)	N104081 (40cm)	ST140-8 3/8 circle (TG140-8)		140°	0.20mm	6.00mm
9-0, 1/2 circle	N96062 (30cm)	N96061 (15cm)	ST160-6 1/2 circle (TG160-6)		160°	0.15mm	5.60mm
9-0, 3/8 circle	N94062 (30cm)	N94061 (15cm)	ST140-6 3/8 circle (TG140-6)		140°	0.15mm	6.00mm
9-0, 3/8 circle	N94082 (40cm)	N94081 (40cm)	ST140-8 3/8 circle (TG140-8)		140°	0.20mm	6.00mm
8-0, 3/8 circle	N84082 (30cm)		ST140-8 3/8 circle (TG140-8)		140°	0.20mm	6.00mm
6-0, 3/8 circle	N640112 (76cm)		ST140-11 3/8 circle		140°	0.28mm	8.00mm

Taper Point





	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED			CURVATURE	WIRE	LENGTH
7-0, 3/8 circle	N74092TP (30cm)*		TP140-9 3/8 circle (BV130-5)		140°	0.23mm	6.00mm

OPHTHALMIC SUTURE


SPECIFICATIONS

POLYGLYCOLIC ACID, VIOLET BRAIDED

Advanced Micro Spatula Needle







	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
6-0, 1/4 circle	A69092 (45cm)*		ST90-9 1/4 circle 	90°	0.23mm	8.0mm
6-0, 1/2 circle	A66092 (45cm)*		ST160-9 1/2 circle 	160°	0.23mm	8.0mm
7-0, 3/8 circle	A74092 (45cm)*		ST140-9 3/8 circle 	140°	0.23mm	6.5mm
8-0, 3/8 circle	A84082 (30cm)*		ST140-8 3/8 circle 	140°	0.2mm	6.0mm

Taper Point

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
8-0, 3/8 circle	A84082TP (30cm)*		TP140-8 3/8 circle 	140°	0.2mm	6.0mm

POLYPROPYLENE, BLUE MONOFILAMENT


Advanced Micro Spatula Needle

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
9-0 straight	P90062 (20cm)		ST000-6 straight 	0°	0.15mm	16.0mm
9-0 , straight and 1/2 circle	P9006062 (30cm)		ST000-6 straight and ST160-6 1/2 circle  	0°, 160°	0.15mm	16.0mm and 6.2m
9-0 straight and 1/4 circle	P9009062 (30cm)		ST000-6 straight and ST90-6 1/4 circle  	0°, 90°	0.16mm	16.0mm
10-0 straight	P100062 (20cm)		ST000-6 straight 	0°	0.15mm	16.0mm

Taper Point

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
10-0, 1/4 circle	P109062TP (20cm)		TP90-6 1/4 circle 	90°	0.15mm	13.0mm

Reverse Cutting



	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
5-0, 3/8 circle	P540151RC (70cm)		RC130-15 3/8 circle 	140°	0.38mm	13.1mm

OPHTHALMIC SUTURE


SPECIFICATIONS

SILK, BLACK BRAIDED


Advanced Micro Spatula Needle

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
8-0, 3/8 circle	S84082 (30cm)	S84081 (15cm)	ST140-8 3/8 circle (TG140-8)	 140°	0.20mm	6.00mm
6-0, 3/8 circle	S640112 (76cm)		ST140-11 3/8 circle	 140°	0.28mm	8.00mm

Taper Point


	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
7-0, 3/8 circle	S74092TP (45cm)*		TP140-9 3/8 circle (BV130-5)	 140°	0.23mm	6.00mm

Cutting


	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
4-0, 3/8 circle		S440171 (76cm)	CT135-17 3/8 circle	 140°	0.43mm	6.00mm

SILK, VIRGIN

Advanced Micro Spatula Needle

	MODEL-DOUBLE ARMED	MODEL-SINGLE ARMED		CURVATURE	WIRE	LENGTH
8-0, 3/8 circle	VS84082 (40cm)	VS84081 (40cm)	ST140-8 3/8 circle (TG140-8)	 140°	0.20mm	6.00mm

* FDA 510(k) only



VISIONARY MEDICAL SUPPLIES, INC.

SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE

6441 Enterprise Lane, Madison, WI 53719 USA TEL: 608-270-3880 FAX: 608-270-3882 www.visionarymedicalsupplies.com