

# VISIFLEXSQ


## Hydrophilic Foldable Intraocular Lenses

VisiFlex lenses are sold with the VisiGlide injector.  
Use VisiGlide to minimize surgical trauma and immediate postoperative inflammation.



NOT AVAILABLE FOR  
SALE IN THE U.S.

### DESCRIPTION



|                   |   |
|-------------------|---|
| Material:         | poly-HEMA                               |
| Refractive Index: | 1.460                                   |
| Optic Diameter:   | 6.0mm                                   |
| Overall Diameter: | 12.0mm                                  |
| Water Content:    | 25%                                     |
| Angulation:       | 0 Degrees                               |
| 'A' Constant*:    | 117.8 to ULIB Standard Measured Post-op |
| 'A' Constant*:    | 118.0 by Manufacturer's Estimate        |
| ACD*:             | 5.0                                     |
| Position Holes:   | 0                                       |

### POWER RANGE

10 to 15 in 1.0 Diopter  
15 to 25 in 0.5 Diopter  
25 to 30 in 1.0 Diopter

### CHARACTERISTICS

- Zero degree angulation facilitates error-free implantation
- Choose incision size: 2.8mm standard; 1.8mm wound assisted
- Faster wound healing and less surgery-induced astigmatism
- Improved unaided acuity and faster visual recovery
- Square edges prevent PCO
- Equiconvex UV absorbing optics for less aberrations
- Clinically proven safety and efficacy
- Specifications comply with ISO 11979 series of standards

\* 'A' Constant and ACD are estimated values only



**VISIONARY MEDICAL SUPPLIES, INC.**  
SUTURES, IOLS, OPHTHALMIC SOLUTIONS, AND MORE

6441 Enterprise Lane, Madison, WI 53719 USA

TEL: 608-270-3880

FAX: 608-270-3882

[www.visionarymedicalsupplies.com](http://www.visionarymedicalsupplies.com)